

4-STAGE ULTRA CLEAN WATER SYSTEM



Quality water Onsite for your Autoclave & Closed Water Systems quickly and affordably

The Ezee Kleen 2.5 HD is designed for offices in geographic areas having tap water mineral levels > 100 PPM (mg/L) or offices using high volumes of water. For geographic areas having tap water mineral levels < 100 PPM (mg/L) or smaller offices consider the Ezee Kleen 2.5 Distillation Alternative.

Contact Oasis Dental Group Inc. for assistance in determining tap water mineral levels and the appropriate unit for your office.

The Ezee Kleen 2.5 HD will:

- Produce 1 gallon (4 litres) of autoclave-safe water in approximately 4 minutes!
- Self monitor for consistent water quality (Unlike bottled water)
- Destroy bacteria, viruses, and other harmful microorganisms
- Reduce suspended solids (organic & inorganic)

- Ensure you never run out of water for your autoclaves
- Eliminate last-minute store runs for water
- Eliminate bulky storage of bottles
- Reduce maintenance and equipment failure
- Save Money!!!
- Produce onsite, water of the quality necessary to comply with all autoclave manufacturer's TDS level guidelines

An Innovative product for the cost effective and efficient dental practice of today.

15531 - 24th Avenue, Unit 24, Surrey, BC V4A 2J3 • Ph: (800) 338-6693 • (604) 538-1285 • Fax: (604) 538-1295

COMPARE EZEEKLEEN 2.5 HD to OnsSide Distillations Alborthed (Distilled) Water	Ezee Kleen 2.5 HD On-Site Distiller		Bottled Distilled Water	
Meets & exceeds autoclave manufacturers specifications	Yes \leq 2.5 PPM (mg/L) ? - possibly, if properly maintained		Inconsistent	
Time required to produce 1 gallon of water	4 minutes	4-8 hours	Time required to pick up water at retail outlet	
Energy requirements Nominal Monitor: 2 x AA Batteries UV Cartridge: 44VAC/6W		4-8 hours of electricity $\approx 5-12 \text{ kWh } (a) 10 \phi = \$0.50-\$1.20/$ gallon ^a	Staff down time - Gas required for vehicle	
Convenience	Convenience Under Counter Installation		Space required to store full & empty bottles	
Microbial Reduction	UV Sterilization at 32,000- 35,000 (w-sec/cm ²)	? - Possibly if properly maintained	Inconsistent	
FEAT	LIRES & S	PECIFICATI	ONS	

WATER ON DEMAND - NEVER RUN OUT! The Ezee Kleen 2.5 HD produces water on demand at a flow rate of 1 gallon (4 litres) in 4 minutes eliminating the need for storage of bottles or runs to the store.

ENERGY EFFICIENT! Monitor powered by 2 x AA Batteries and a 44VAC/6W UV Cartridge transformer the Ezee Kleen 2.5 HD is extremely energy efficient operating for pennies a day. Generates no heat to present hazardous fire concerns or discomfort to staff and patients.

MAXIMIZE EFFICIENCY AND LIFE OF AUTOCLAVE EQUIP-MENT! The Ezee Kleen 2.5 HD comes complete with 1 Prefiltration cartridge, 1 RO cartridge, 1 mixed bed DI resin cartridge and 1 UV cartridge that bring your TDS levels to as low as 0.4 PPM (mg/L) ensuring your equipment does not become clogged with mineral scale that can cause malfunction, increased energy usage and shortened life. (Autoclave manufacturers recommend water quality of less than 5.0 PPM (mg/L) is ideal for maximum efficiency and life.)

PROVIDE PUREST OF WATER POSSIBLE FOR CONTAINED DEN-TAL UNIT WATER SYSTEMS!

No microbes means less chance of biofilm growth to begin in the first place, slower growth if it does occur. Reduce chemical treatment frequency on your dental lines to save time, money and deterioration of lines and equipment from exposure to harsh chemicals. Reduce the change frequency of micro-filtration cartridges on your dental lines saving time and money.

FAST SIMPLE INSTALLATION AND USE! Unit comes ready for use with all necessary components for under counter installation. Estimated installation time: 2 hours. **SECONDS TO CHANGE CAR-TRIDGES & NO MESS!** Staff members can easily change cartridges. Cartridges are self contained sealed units with quick connect fittings allowing easy cartridge change out that take seconds and no exposure to contaminants.

APPROVED COMPONENTRY AND ONE YEAR WARRANTY! All components are UL/CSA approved, made of FDA approved material, and require no MSDS reporting. The Ezee Kleen 2.5 HD comes complete with one year parts and labour (excluding cartridges) return to depot warranty.

NO MORE GUESSING WATER QUALITY! The Ezee Kleen 2.5 HD self monitoring meter presents a digital read out of water quality indicating when DI Cartridge change is required.

Gallon Production ^b of EZEEKLEEN 2.5 HD	DI Cartridge
based on PPM* (mg/L) of TDS** in	Tap Water

TDS Levels	Production US Gallons (Litres)			Production with XL DI US Gallons (Litres)	
100 (1)	529	(2000)	2192	(8285)	
200	264	(1000)	1095	(4143)	
400 (2)	132	(500)	548	(2071)	
800	66	(250)	274	(1036)	

* PPM = Parts Per Million = mg/L **TDS = Total Dissolved Solids Actual gallon production from EZEEKLEEN 2.5 HD DI Cartridge will depend on TDS levels in source tap water and other contaminants.

levels in source tap water and other contaminants. (1) Ezee Kleen 2.5 Distillation Alternative is more cost efficient in areas where TDS < 100 PPM (mg/L)

(2) In geographic areas where TDS is > 400 PPM (mg/L) use of the 8002-XL DI Cartridge will reduce DI Cartridge change frequency and improve cost per gallon.

Unit Options

8000 Ezee Kleen 2.5 HD - 2 x AA Batteries & 44VAC/ 6W North American Format

8000-XL Ezee Kleen 2.5 HD XL - for geographic areas having tap water dissolved mineral levels > 400 mg/L; Complete with 8002-XL-KIT

8000-14 Ezee Kleen 2.5 HD with 14 Gallon Tank - for large offices going through high volumes of water. (Contact Oasis Dental Group Inc. for tank sizes in excess of 14 gallons.)

8000+RO Ezee Kleen 2.5 HD with RO System - Provides Stage 2 - RO Water for the rinse cycle of instrument washers, SCWS's, etc.

Replacement Cartridges

8002-FR DI Cartridge - mixed bed DI resins remove dissolved solids and bring the quality of the water to as low as .04 PPM of total dissolved solids. Change when monitor indicates (> $5.0 \ \mu$ S with water flowing.)

8003 Prefiltration Cartridge - granular activated carbon filter which removes impurities that affect the aesthetic quality of water and prepares the water for RO processing. Change annually.

8004 RO Cartridge - removes 90-95% of organic material such as algae, bacteria, amoebae etc. and suspended solids such as sand, silt, metal particles (rust flakes) etc.. Change bi-annually (every two years)

8005 UV Cartridge - High powered , energy efficient 32,000 w-sec/cm² ultraviolet radiation destroys bacteria, viruses and other micro-organisms by interfering with the DNA and RNA in the organisms reproductive cycle. Change annually.

Cartridge Options

Recommended for geographic areas having tap water dissolved mineral levels > 400 mg/L or for offices using high volumes of water. Change with 8002-FR when monitor indicates (> 5.0 μ S with water flowing.)

8002-XL-KIT Extra Large DI Cartridge Kit - Contains all parts to install XL DI Cartridge for the first time (cartridge included.) Replace with 8002-XL

8002-XL Extra Large DI Cartridge - Reduce change frequency of Ezee Kleen 2.5 HD DI Cartridge (8002-FR) by as much as 5 times and improve cost/gal



Unit Dimensions: 16"H x 12"W x 6"D/ 41cm x 30cm x 15cm Tank Dimensions: 12" D x 17" H/ 31cm D x 43cm H Weight: 27 lbs/ 12 kg

a. Information is for comparative purposes only and believed to be accurate based on calculations by Oasis Dental Group Inc., July 2011. Actual energy costs will vary. b. These data are believed to be reliable but do not imply any warranty or performance guarantee. Actual performance will vary. Microbial testing may be appropriate to determine suitability of water for intended use.