

# EZEKLEEN 2.5 HD

4 - STAGE ULTRA CLEAN WATER SYSTEM  
**Optional XL DI Cartridge**



## Increase the Life of Your Ezee Kleen 2.5 HD DI Cartridge by 4-5 Times!

*Designed for offices wanting to reduce the change frequency of the 8002-FR Ezee Kleen 2.5 HD DI Cartridge*

### ***Change Ezee Kleen 2.5 HD DI Cartridge less frequently***

The Ezee Kleen 2.5 HD XL DI Cartridge packs tons of Deionization resin extending the life of the Ezee Kleen 2.5 HD DI Cartridge by 4-5 times!

### ***Spend Less Money on Ezee Kleen 2.5 HD DI Cartridges***

- The Ezee Kleen 2.5 HD XL DI Cartridge packs 4-5 times more Ezee Kleen 2.5 HD DI resin into a larger housing at a significantly reduced price!

### ***Save Service Costs & Time***

### ***Ideal For:***

- Areas having high TDS levels
- Offices using high volumes of autoclave water

### ***The Ezee Kleen 2.5 HD XL DI Cartridge***

- Order as an option with a new Ezee Kleen 2.5 HD
- Easily retrofitable into existing Ezee Kleen 2.5 HD's
- Change when monitor indicates along with the 8002-FR DI Cartridge

### ***Fast Simple Installation!***

- Take minutes to install
- No special tools required

**Ezee Kleen 2.5 HD XL DI Install Kit:** 8002-XL-KIT (Includes one XL DI Cartridge)

**Ezee Kleen 2.5 HD XL DI Replacement Cartridge:** 8002-XL (Replaced when monitor Indicates)



15531 - 24th Avenue, Unit 24, Surrey, BC V4A 2J3 • Ph: (800) 338-6693 • (604) 538-1285 • Fax: (604) 538-1295

# EZEKLEEN 2.5 HD

4 - STAGE ULTRA CLEAN WATER SYSTEM

Model: 8000

## 8002-XL DI CARTRIDGE INSTALL DIAGRAM


### Dimensions:

Unit - Length 12" x Height 16" x Depth 6"

Tank - Diameter 11" x Height 17"

8002-XL DI Cartridge + Receiving Cap + Mounting Bracket - Base Diameter 3" x Head Diameter 4" x Height 20"

XL DI MOUNTING BRACKET →



**Oasis**  
DENTAL GROUP INC.

15531 - 24th Avenue, Unit 24, Surrey, BC V4A 2J3 • Ph: (800) 338-6693 • (604) 538-1285 • Fax: (604) 538-1295