

MAINTENANCE SUPPLIES

TO ORDER YOUR SUPPLIES

1) Contact your authorized **EZEEKLEEN 2.5** dealer.

For easy reference, enter your Dealer's name and contact information here:

Or

2) Contact **Oasis Dental Group** at:

1-800-338-6693 or 1-604-538-1285

Fax: 1-604-538-1295

Website: www.oasisdentalgroup.ca

Email: info@oasisdentalgroup.ca

EzeeKleen 2.5 - Complete Unit - # 9000

De-ionization Cartridges - # 9001

- change cartridges when monitor reads greater than 5.0 μS **with water flowing** water is greater than 2.5 PPM (mg/L) of dissolved solids and is unsafe to use in autoclave.

OPTIONAL FEATURES

- **XL DI Cartridge - #9001-XL**
Change at the same time as 9001 DI Cartridges (2 Pak) when monitor reads greater than 5.0 μS with water flowing
- **XL DI Cartridge Install Kit - 9001-XL-KIT**
Order for initial XL DI Cartridge install.
- **Sediment Cartridge - #9002**
Change after one year of use.
- **Sediment Cartridge Install Kit - #9003**
Order for initial Sediment Cartridge install.
- **Ultra-Violet Cartridge - #9004**
Order for initial Ultra-Violet Cartridge install.
- **Ultra-Violet Cartridge Replacement Lamp - #9004-1**
Change after one year of use.
- **Under Sink Installation Kit - #9005** (see page. 8)
- **Permanent Mount Faucet - #9007** (see page. 9)
- **Retractable Faucet Kit - #8020**
- Total Dissolved Solids Tester (Handheld)
- For water testing contact Oasis Dental Group or a local water testing facility.

TIP: Use Patient scheduler to log a reminder to change Sediment Cartridge and Ultra Violet Lamp after one year of use.



15531 - 24th Avenue, Unit 24, Surrey, BC V4A 2J3 • Ph: (800) 338-6693 • (604) 538-1285 • Fax: (604) 538 - 1295